



Recombinant Bacillus pseudofirmus Na (+)/H (+) antiporter subunit E (mrpE)

Product Code	CSB-CF865687BRB
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9RGZ1
Product Type	Transmembrane Protein
Immunogen Species	Bacillus pseudofirmus (strain OF4)
Sequence	MAFQILLNLVIAVIWVNFQNSYTAVDFLIGYVVGIFILFVLRRLRFLRDFYMRRIWAI IKL ISLFFKELILANIDVIKIVLSPKMNIQPGIVAVPTKLKTDWELSLLASLISLTPGTLS MD FSDDNKYIYIHAIIDVPNKEK MIRDIDHDTFERAILEVTK
Source	in vitro E.coli expression system
Target Names	mrpE
Protein Names	Recommended name: Na(+)/H(+) antiporter subunit E Alternative name(s): Mrp complex subunit E Multiple resistance and pH homeostasis protein E
Expression Region	1-158
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.