



Recombinant UPF0041 protein R07E5.13 (R07E5.13)

Product Code	CSB-CF638610CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q21828
Product Type	Transmembrane Protein
Immunogen Species	Caenorhabditis elegans
Sequence	MSRVISKVTTYFKQHSTAEWKHYFLSTHFWGPVANWGLPLAALGDLKKNPDM ISGPMTSA LLIYSSVFMRFQVHWQPRNLLLFACHFANFSAQGAQLGRFVNHNYLHYVEDPV HHKLMMK KEVLEHEHDAERVSLFLWYPFFEFLGNFRKNKLNKNKFM
Source	in vitro E.coli expression system
Target Names	R07E5.13
Protein Names	Recommended name: UPF0041 protein R07E5.13
Expression Region	1-160
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.