



Recombinant Danio rerio Protein cornichon homolog 2 (cnih2)

Product Code	CSB-CF710680DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q5BL21
Product Type	Transmembrane Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Sequence	MAFTFAAFCYMLTLVLC AALIFFVIWQIIAFDELRTDFKNPIDQSNPTRARERILNI ERI CNLLRRLVVPEYSIHGLFCLMFMCAGEWVTLGLNIPLLLYHLWRFFHRPADGS EVMYDPV SVMNADILNYCQKESWCKLGFYLLSFFYYLYSMVYALVSF
Source	in vitro E.coli expression system
Target Names	cnih2
Protein Names	Recommended name: Protein cornichon homolog 2
Expression Region	1-160
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.