



Recombinant Chlamydomonas reinhardtii Sulfur-rich protein (srp)

Product Code	CSB-CF310171DSL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P94665
Product Type	Transmembrane Protein
Immunogen Species	Chlamydomonas reinhardtii (strain GPIC)
Sequence	MLTGVENSESGVIDLIKPLDDVMKNETVQVTLVNSVLGWCKAHIVDPIKTSKI VQSRAF QITMVVLGVILLIAGLALTFVLQGGQLGKNAFLFLIPAVIGLVKLLTTSVFMEKPCTP EKW RLCKRLLATTEDILDDGQINQSNTIFTTESSDVTNTATQS
Source	in vitro E.coli expression system
Target Names	srp
Protein Names	Recommended name: Sulfur-rich protein
Expression Region	1-160
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.