



Recombinant African swine fever virus Protein MGF 110-13L (War-020)

Product Code	CSB-CF315302AEI
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0C9K0
Product Type	Transmembrane Protein
Immunogen Species	African swine fever virus (isolate Warthog/Namibia/Wart80/1980) (ASFV)
Sequence	MGGGDYWPIIIRHCCFYLVFSIAFVG YIVFAYYKNLHLNNTMKLIALLCILIWLSQ PGLN RPLSIFYMKQNLPRTYTPPIRELEYWCTYGKHCDFCWECRNGICKKNKVWDDM PLIKQNDY ISQCSIARYFDRCMYFIKPKTPYIHYMDCSQPTAYKGF SH
Source	in vitro E.coli expression system
Target Names	War-020
Protein Names	Recommended name: Protein MGF 110-13L
Expression Region	1-160
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.