



Recombinant *Saccharomyces cerevisiae* UPF0479 membrane protein YER190C-B (YER190C-B)

Product Code	CSB-CF317849SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0CX94
Product Type	Transmembrane Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Sequence	MMPAKLQLDVLRTLQSSARHGTQTLKNSNFLERFHKDRIVFCLPFFPALFFVPV QKVLQH LCLRFTQVAPYFIIQLFDLPSRHAENLAPLLASCRIQYTNCFSSSSNGQVPSIISL YLRV DLSPFYAKIFQISYRVPMIWLDFVQVFFVFLVISQHSLS
Source	in vitro E.coli expression system
Target Names	YER190C-B
Protein Names	Recommended name: UPF0479 membrane protein YER190C-B
Expression Region	1-160
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.