



Recombinant *Saccharomyces cerevisiae* Protein SLM4 (SLM4)

Product Code	CSB-CF336424SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P38247
Product Type	Transmembrane Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Sequence	MVMLHSKNVKGFLENTLKPVDLHSDVDFKTSSLQSSMIITATNGGILSYATSNND VPKNSI NEINSVNNLKMMSLLIKDKWSEEDNDTEEQHSNSCYPVEIDSFKTKIYTYEMED LHTCVA QIPNSDLLLFFIAEGSFPYGLLVKIERAMRELDFGYKLG
Source	in vitro E.coli expression system
Target Names	SLM4
Protein Names	Recommended name: Protein SLM4 Alternative name(s): EGO complex subunit 3 GSE complex subunit 1
Expression Region	1-162
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.