



# Recombinant *Macaca fascicularis* NADH dehydrogenase [ubiquinone] 1 alpha subcomplex subunit 13 (NDUFA13)

<b>Product Code</b>	CSB-CF700164MOV
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	Q4R6H1
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	<i>Macaca fascicularis</i> (Crab-eating macaque) ( <i>Cynomolgus</i> monkey)
<b>Sequence</b>	MAVAVCHFRLGPEVWNTASMEMPKVKQDMPPGGYGPIDYKRNLPRRGLSG YSMLAIGIG TLVYGHWSIMKWNRRERRRLQIEDFEARIALMPLFQAETDRRTLQMLRENLEEE AIIMKDV PDWKVGESVFHTTRWVPLIGELYGLRTTEETIHANYGFMWYT
<b>Source</b>	in vitro <i>E.coli</i> expression system
<b>Target Names</b>	NDUFA13
<b>Protein Names</b>	Recommended name: NADH dehydrogenase [ubiquinone] 1 alpha subcomplex subunit 13 Alternative name(s): Complex I-B16.6 Short name= CI-B16.6 NADH-ubiquinone oxidoreductase B16.6 subunit
<b>Expression Region</b>	1-163
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.