



Recombinant *Saccharomyces cerevisiae* Protein NAG1 (NAG1)

Product Code	CSB-CF855084SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8TGN9
Product Type	Transmembrane Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Sequence	MNSAGRVHRSRAGSRGHAAISPLTMASFSVARGIRSSNVYDDTDDELSILTFF SAVRRNR LTSSLPPILSARCSSACFSVRIVLPLSLTISISALMYSTNSALGRKLTGAFSIQTNI EQS CGFFRTSIMATLPPIECPPIIGPPLVFNSCFVIKCFSSDMS
Source	in vitro E.coli expression system
Target Names	NAG1
Protein Names	Recommended name: Protein NAG1 Alternative name(s): Nested antisense gene 1 protein
Expression Region	1-163
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.