



Recombinant Escherichia coli O139:H28 UPF0114 protein YqhA (yqhA)

Product Code	CSB-CF426105EJD
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A7ZRN9
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli O139:H28 (strain E24377A / ETEC)
Sequence	MERFLENAMYASRWLLAPVYFGLSLALVALALKFFQEIIHVLPNIFSLMAESDLIL VLLSL VDMTLVGGLLVMVMFSGYENFVSQLDISENKEKLNWLGKMDATSLKNKVAASI VAISSIH LLRVFMDAKNVPDNKLMWYVIIHLTFVLSAFVMGYLDRLTRHNH
Source	in vitro E.coli expression system
Target Names	yqhA
Protein Names	Recommended name: UPF0114 protein YqhA
Expression Region	1-164
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.