



Recombinant Human Glycosylphosphatidylinositol-anchored high density lipoprotein-binding protein 1 (GPIHBP1)

Product Code	CSB-CF811603HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8IV16
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	QTQQEEEEEDHGPDDYDEEDEDEVEEEETNRLPGGRSRVLLRCYTCKSLP RDERCNLT QNC SHGQTCTTLIAHGNTESGLLTTHSTWCTDSCQPITKTVEGTQVTMTCCQ SSLCNVPP WQSSRVQDPTGKGAGGPRGSSETVGAALLLNLLAGLGAMGARRP
Source	in vitro E.coli expression system
Target Names	GPIHBP1
Protein Names	Recommended name: Glycosylphosphatidylinositol-anchored high density lipoprotein-binding protein 1 Short name= GPI-HBP1 Short name= GPI-anchored HDL-binding protein 1 Alternative name(s): High density lipoprotein-binding protein 1
Expression Region	21-184
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.