



Recombinant Human Uncharacterized protein C3orf45 (C3orf45)

Product Code	CSB-CF822677HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8N112
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MPSLAPDCPLLAMPEETQEDSVAPMMPSQRSRGPLAPNHVHEVCLHQVESIS DLHSGAGT LRPYLTEEARPWDELLGVLPPSLCAQAGCSPVYRRGGFLLLLALLVLTCLVLAL LAVYLS VLQSESLRILAHTLRTQEETLLKLRLASLSQLRRLNSSEAQAPS
Source	in vitro E.coli expression system
Target Names	LSMEM2
Protein Names	Recommended name: Uncharacterized protein C3orf45
Expression Region	1-164
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.