



Recombinant Danio rerio Mitochondrial fission process protein 1 (mtfp1)

Product Code	CSB-CF764780DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q6PCS6
Product Type	Transmembrane Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Sequence	MEPSADTHSKPVDIYRDTWVRFLGYANEVGEAFRALVPVGAVWASYGVATTY VTADAIDK GRKAAAAHGERPGKAVCVCAVVDTFVWQALASVAVPGFTINRVCAASHFLL SRTTRWPL PVRKWTTTAIGLSTIPFIITPIDRSVDLLLDSSLRKLYSEGEKED
Source	in vitro E.coli expression system
Target Names	mtfp1
Protein Names	Recommended name: Mitochondrial fission process protein 1 Alternative name(s): Mitochondrial 18 kDa protein Short name= MTP18
Expression Region	1-165
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.