



Recombinant Emerin homolog 1 (emr-1)

Product Code	CSB-CF519041CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O01971
Product Type	Transmembrane Protein
Immunogen Species	Caenorhabditis elegans
Sequence	MDVSQLTDAELRDSLKSHGVS VGPVATTRKLYEKKLIKLS DGSINNQSNL NDS QFNEDS LISSSPKKSPQRVFQNVSAATAAATTSPESDSDDC EESMRYL TEEEMAADR ASARQAQ SNKGGFLGSTITFTILFVFI AVFAYFLIENAEQLKLV AETNPEDI
Source	in vitro E.coli expression system
Target Names	emr-1
Protein Names	Recommended name: Emerin homolog 1 Alternative name(s): Ce-emerin
Expression Region	1-166
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.