



Recombinant Zea mays Cell number regulator 3 (CNR3)

Product Code	CSB-CF515173ZAX
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	D9HP19
Product Type	Transmembrane Protein
Immunogen Species	Zea mays (Maize)
Sequence	MYPATTPYETASGVGVAPVAGLFPVAGEAREWSSRLLDCFDDFDICCMTFWC PCITFGRT AEIVDHGMTSCG TSAALFALIQLWLSGSQCTWAFSCTYRTRLRAQHGLPEAPC ADFLVHLC CLHCALCQEYRELKARGYEPVLGWEFNAQRAAAGVAMCPPASQGMGR
Source	in vitro E.coli expression system
Target Names	CNR3
Protein Names	Recommended name: Cell number regulator 3 Alternative name(s): ZmCNR03
Expression Region	1-167
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.