



Recombinant *Ashbya gossypii* Genetic interactor of prohibitin 7, mitochondrial (GEP7)

Product Code	CSB-CF758786DOT
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q75EB9
Product Type	Transmembrane Protein
Immunogen Species	<i>Ashbya gossypii</i> (strain ATCC 10895 / CBS 109.51 / FGSC 9923 / NRRL Y-1056) (Yeast) (<i>Eremothecium gossypii</i>)
Sequence	SDTVSKLNSGELTEVVRQRYCVDNKRCEMRMLLTKYPGPAREREMQLQVATA ELSDWR KMPRVWQQVSYHAFGSWGPRTGLSFVGRRPEDFFVTDQKGLWTCSAPRR AEYERSSRAL DPASRAVLYAAALVAAVAALGDLWRRQDADRQVTVAELDLAEPSSPT
Source	in vitro E.coli expression system
Target Names	GEP7
Protein Names	Recommended name: Genetic interactor of prohibitin 7, mitochondrial
Expression Region	67-234
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.