



Recombinant *Saccharomyces cerevisiae* Putative uncharacterized transporter YOL163W (YOL163W)

Product Code	CSB-CF316562SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0CF20
Product Type	Transmembrane Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Sequence	MIAWSLVATLQCKMTGKSSFYTCRALMGLFEGGFVADLVLMWSYFYSSSELSI RLSFFWV TSLTQIITSIVAFGVFHMARGIGGMAGWQWLFLIERIFTLVIGISAYFLMVPSVVQ TKKP WSKKGWFTEREEKIIVNKILRDDPTKGDMMNRQGMSLKMLWQGITDYYI
Source	in vitro <i>E.coli</i> expression system
Target Names	YOL163W
Protein Names	Recommended name: Putative uncharacterized transporter YOL163W
Expression Region	1-169
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.