



Recombinant Clostridium botulinum UPF0316 protein CLJ_B0679 (CLJ_B0679)

Product Code	CSB-CF501389DUE
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	C3L101
Product Type	Transmembrane Protein
Immunogen Species	Clostridium botulinum (strain 657 / Type Ba4)
Sequence	MLSYYAFFFFAKIMEVALMTIRTVLITRGEKLYGSIIGFIEVTIWLYVTSSVLSGIKD DP IRMVVYALGFTCGNYMGCVIEEKLAIGLLTINVITSESDGKRLAEILRDKNVGVT MVDAE GKIEQKKMLIIHAKRKRREEIIRTIEGSDINAMISVNDIKTVYGGYGIRK
Source	in vitro E.coli expression system
Target Names	CLJ_B0679
Protein Names	Recommended name: UPF0316 protein CLJ_B0679
Expression Region	1-170
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.