



Recombinant Human TM2 domain-containing protein 1 (TM2D1)

Product Code	CSB-CF887145HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9BX74
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	TSAGGEESLKCEDLKVGQYICKDPKINDATQEPVNCTNYTAHVSCFPAPNITCK DSSGNE THFTGNEVGFFKPISCRNVNGYSYKVAVALSLFLGWLGADRFYLGYPALGLLK FCTVGFC GIGSLIDFILISMQIVGPSDGSSYIIDYYGTRLTRLSITNETFRKTQLYP
Source	in vitro E.coli expression system
Target Names	TM2D1
Protein Names	Recommended name: TM2 domain-containing protein 1 Alternative name(s): Beta-amyloid-binding protein Short name= hBBP
Expression Region	38-207
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.