



Recombinant *Oryza sativa* subsp. *indica* Oleosin 18 kDa (OLE18)

Product Code	CSB-CF386080OFF
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A2XL05
Product Type	Transmembrane Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>indica</i> (Rice)
Sequence	ADRDRAGQYYQQQRGQVGETVKGILPEKAPSASQALTVATLFPGLGLLLLVLSG LALAASV VGLAVATPVFLIFSPVLVPAALLIGLAVAGFLTSGALGLGGLSSLTFLANTARQA FQRTF DYVEQARRRMAEAAAHAAGHKTAQAGHAIQGRADQAGTGAGAGGGAGTKTS S
Source	in vitro E.coli expression system
Target Names	OLE18
Protein Names	Recommended name: Oleosin 18 kDa Alternative name(s): OSE721
Expression Region	2-172
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.