



Recombinant *Drosophila melanogaster* Transmembrane protein 70 homolog, mitochondrial (CG7506)

Product Code	CSB-CF856836DLU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q95SS8
Product Type	Transmembrane Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Sequence	SDDALQRIYYGTLAPRMKMKVFFSLSTSLAGLAAQPILLEQGMKIGGTGMAVFL CTVGGF FTFVTPLLLHFITKKYVTELHYNPLTEEYATTISLLLQKIKTTFRPNDVVVPEVP GMFT SFLVNRPLFVDPALFDDPEHYVKIMGYDKPIDFKLDLTPNPEKSKTSEEKQ
Source	in vitro <i>E.coli</i> expression system
Target Names	CG7506
Protein Names	Recommended name: Transmembrane protein 70 homolog, mitochondrial
Expression Region	65-236
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.