



Recombinant Human Calcium channel flower homolog (CACFD1)

Product Code	CSB-CF880082HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9UGQ2
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MSSSSGGAPGASASSAPPAQEEGMTWWYRWLCRLSGVLGAVSCAISGLFNCI TIHPLNIAA GVWMIMNAFILLCEAPFCCQFIEFANTVAEKVDRLRSWQKAVFYCGMAVVPI VISLTLT TLLGNAIAFATGVLYGLSALGKKGDAISYARIQQQRQQADEEKLAETLEGEL
Source	in vitro E.coli expression system
Target Names	CACFD1
Protein Names	Recommended name: Calcium channel flower homolog Alternative name(s): Calcium channel flower domain-containing protein 1
Expression Region	1-172
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.