



Recombinant Human Mitochondrial import inner membrane translocase subunit Tim17-B (TIMM17B)

Product Code	CSB-CF023551HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O60830
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MEEYAREPCPWRIVDDCGGAFTMGVIGGGVFQAIKGFRNAPVGIRHRLRGSA NAVRIRAP QIGGSFAVWGGLFSTIDCGLVRLRGKEDPWNSITSGALTGAVLAARSGPLAMV GSAMMGG ILLALIEGVGILLTRYTAQQFRNAPPFLEDPSQLPPKDGTAPGYPSYQQYH
Source	in vitro E.coli expression system
Target Names	TIMM17B
Protein Names	Recommended name: Mitochondrial import inner membrane translocase subunit Tim17-B
Expression Region	1-172
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.