



Recombinant Escherichia coli Phosphatidylglycerophosphatase A (pgpA)

Product Code	CSB-CF324695ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P18200
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli (strain K12)
Sequence	MTILPRHKDVAKSRLKMSNPWHLLAVGFGSGLSPIVPGTMGSLAAIPFWYLMT FLPWQLY SLVVMLGICIGVYLCHQTAKDMGVHDHGSIVWDEFIWMWITLMALPTNDWQW VAAGFVIF RILDMWKWPPIRWFDRNVHGGMGIMIDDIVAGVISAGILYFIGHHWPLGILS
Source	in vitro E.coli expression system
Target Names	pgpA
Protein Names	Recommended name: Phosphatidylglycerophosphatase A EC= 3.1.3.27
Expression Region	1-172
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.