



Recombinant Human Tail-anchored protein insertion receptor WRB (WRB)

Product Code	CSB-CF026148HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O00258
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MSSAAADHWAWLLVLSFVFGCNVLRILLPSFSSFMSRVLQKDAEQESQMRAEI QDMKQEL STVNMMDEFARYARLERKINKMTDKLKTHVKARTAQLAKIKWVISVAFYVLQAA LMISLI WKYYSPVAVVPSKWITPLDRLVAFPTRVAGGVGITCWILVCNKVVAIVLHPFS
Source	in vitro E.coli expression system
Target Names	WRB
Protein Names	Recommended name: Tail-anchored protein insertion receptor WRB Alternative name(s): Congenital heart disease 5 protein Tryptophan-rich basic protein Short name= WRB
Expression Region	1-174
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.