



Recombinant TM2 domain-containing protein C41D11.9 (C41D11.9)

Product Code	CSB-CF836118CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q95QZ5
Product Type	Transmembrane Protein
Immunogen Species	Caenorhabditis elegans
Sequence	ISGTKDVKSKNCDGSAGLTCTFPGD CRIGDTVKNCTSRKGCNPVSRNNVE AVCRFCWQ LLPGDYDCEPATNCSTSSTKLLVTKCSAHSSVICMGQRNFYKRIPCNWSSGYS WTKTML SVVLGGFGADRFYLGLWKS AIGKLFSFGGLGVWTLVDVVLIAVGYIKPYDGSM YI
Source	in vitro E.coli expression system
Target Names	C41D11.9
Protein Names	Recommended name: TM2 domain-containing protein C41D11.9
Expression Region	21-195
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.