



Recombinant Human Tumor suppressor candidate 5 (TUSC5)

Product Code	CSB-CF025354HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8IXB3
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MAHPVQSEFPQAQEPGSAAFDLPEMEILLTKAENKDDKTLNLSKTLGSLDLE QNSQGL PFKAISEGHLEAPLPRSPSRASSRRASSIATTSYAQDQEAPRDYLLAVVACFC PVWPLN LIPLIISIMSRSSMQQGNVDGARRLGRLARLLSITLIIMGIVIIMVAVTVNFTVQKK
Source	in vitro E.coli expression system
Target Names	TUSC5
Protein Names	Recommended name: Tumor suppressor candidate 5 Alternative name(s): Interferon-induced transmembrane domain-containing protein D3 Protein located at seventeen-p-thirteen point three 1
Expression Region	1-177
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.