



Recombinant Human Transmembrane protein 52 (TMEM52)

Product Code	CSB-CF023852HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8NDY8
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	FADGSCDPSDQCPPQARWSSLWHVGLILLAVLLLLLLCGVTAGCVRFCCLRKQ AQAQPHLP PARQPCDVAVIPMDSDSPVHSTVTSYSSVQYPLGMRLPLPFGELDLDSMAPP AYSLYTPE PPPSYDEAVKMAKPREEGPALSQKPSPLLGASGLETTVPQESGPNTQLPPC SPGAP
Source	in vitro E.coli expression system
Target Names	TMEM52
Protein Names	Recommended name: Transmembrane protein 52
Expression Region	33-209
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.