



Recombinant Chicken Anti-apoptotic protein NR13 (NR13)

Product Code	CSB-CF852696CH
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q90ZN1
Product Type	Transmembrane Protein
Immunogen Species	Gallus gallus (Chicken)
Sequence	MPGSLKEETALLLEDYFQHRAGGAALPPSATAAELRRAAAELERRERPFFRSC APLARAE PREAAALLRKVAAQLEAEGGLNWGRLLALVVFTGTLAAALAESGCEEGPSRLA AALAAYL AEEQGEWLEEHGGWDGFCRFFGRHGSQPADQNSTLSNAIMAAAGFGIAGLA FLLVVR
Source	in vitro E.coli expression system
Target Names	NR13
Protein Names	Recommended name: Anti-apoptotic protein NR13 Alternative name(s): Apoptosis regulator Nr-13
Expression Region	1-177
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.