



Recombinant Arabidopsis thaliana HVA22-like protein a (HVA22A)

Product Code	CSB-CF882792DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9S7V4
Product Type	Transmembrane Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Sequence	MGSGAGNFLKVLLRNFDVLAGPVVSLVYPLYASVQAIETQSHADDKQWLTYW VLYSLLTL IELTFAKLI EWLP IWSYMKLILTCWLVIPYFSGAAYVYEHFVRP VFVNPR SINIWY V P K K MDIFRKPDDVLTAAEKYIAENGPDAFEKILSRADKSKRYNKHEYESYETMYGE GYQY
Source	in vitro E.coli expression system
Target Names	HVA22A
Protein Names	Recommended name: HVA22-like protein a Short name= AtHVA22a
Expression Region	1-177
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.