



Recombinant Danio rerio Lipoma HMGIC fusion partner homolog (lhfp)

Product Code	CSB-CF732512DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q5PRC1
Product Type	Transmembrane Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Sequence	TSCVGGFFMPYWLLGSQMGPVSFGTFRRCSEYPIRDEARGGTVMLEQCGRYA SFQGIPSLE WRICTVVTGIGCGLLLLLVALTAIMGCCVTDLISRTIGRVAGGIQFVGGLLIGSGC ALYPL GWDSEEVQRQTCSNSSDQFDLGSCEIGWAYYYCTGAGAAAAMVLCTWMACTA GKKQKHYPY
Source	in vitro E.coli expression system
Target Names	lhfp16
Protein Names	Recommended name: Lipoma HMGIC fusion partner homolog
Expression Region	22-200
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.