



Recombinant Mouse E3 ubiquitin-protein ligase RNF5 (Rnf5)

Product Code	CSB-CF019894MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O35445
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	AAAEEDGGPEGPNRERGGASATFECNICLETAREAVVSVCGHLYCWPCCLHQ WLETRPDR QECPVCKAGISREKVVPLYGRGSQKQPRLKTPPRPQGQRPAPESRGGFQ PFGDAGGFH FSFGVGAFPFGFFTTVFNAHEPFRRGAGVDLGQGHPASSWQDSLFLFLAIFFF FWLLSI
Source	in vitro E.coli expression system
Target Names	Rnf5
Protein Names	Recommended name: E3 ubiquitin-protein ligase RNF5 EC= 6.3.2.- Alternative name(s): RING finger protein 5
Expression Region	2-180
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.