



Recombinant Human LEM domain-containing protein 1 (LEMD1)

Product Code	CSB-CF715035HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q68G75
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MVDVKCLSDCKLQNQLEKLGFSPPILPSTRKLYEKLVQLLVSPPCAPPVMN GPRELDG AQDSDDSEELNIIHQNIILSTEKSKKLLKKWPEASTTKRKAVDYCLDYKPSKG RRWAAR APSTRITYGTITKERDYCAEDQTIESWREEGFPVGLKLAVLGIFIIVVFVYLTVEN KSLF G
Source	in vitro E.coli expression system
Target Names	LEMD1
Protein Names	Recommended name: LEM domain-containing protein 1 Alternative name(s): Cancer/testis antigen 50 Short name= CT50 LEM domain protein 1 Short name= LEMP-1
Expression Region	1-181
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.