



Recombinant Bacillus anthracis UPF0316 protein BAA_3454 (BAA_3454)

Product Code	CSB-CF503313BQF
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	C3P314
Product Type	Transmembrane Protein
Immunogen Species	Bacillus anthracis (strain A0248)
Sequence	MLQALLIFVLQIIYVPILTIRTILLVKNQTRSAAAVGLLEGAIYIVSLGIVFQDLSNW MN IVAYVIGFSAGLLLGGYIENKLAIGYITYQVSLLDRCNELVDELRHSGFGVTVFE GEGIN SIRYRLDIVAKRSREKELLEIINEIAPKAFMSSYEIRSFKGGYLTAMKKRALMKK KDHH VS
Source	in vitro E.coli expression system
Target Names	BAA_3454
Protein Names	Recommended name: UPF0316 protein BAA_3454
Expression Region	1-182
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.