



Recombinant Arabidopsis thaliana Oleosin 5 (At3g01570)

Product Code	CSB-CF871146DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9SS98
Product Type	Transmembrane Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Sequence	ADVRTHSHQLQVHPQRQHEGGIKVLYPQSGPSSTQVLAVFVGVPPIGGTLLTIA GLTLAGS VIGLMLAFPLFLIFSPVIVPAAFVIGLAMTGFLASGAIGLTGLSSMSWVLNYIRRA GQHI PEELEEAHRLADMAEYVGQRTKDAGQTIEDKAHDVREAKTFDVRDRDRTTKG THNVRDTK TT
Source	in vitro E.coli expression system
Target Names	At3g01570
Protein Names	Recommended name: Oleosin 5
Expression Region	2-183
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.