



Recombinant Uncharacterized protein F46C5.2 (F46C5.2)

Product Code	CSB-CF346026CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P52881
Product Type	Transmembrane Protein
Immunogen Species	Caenorhabditis elegans
Sequence	QMTCEEFELENTESLVECVQCVQLWRSASAKEGKVCTSGATTCKGNACFMR QCKHCPVYQ YMSGCVNFSPWQLADLEMNRRRTSELRMRRVGAVLLCEDTFNQTTVCNRRD KCNDIHSRL PFATYAEGFLFRGVVNFDTIIAAIDPRYLEVMSGYHFRFLASSSSSSFSFLPSIAILL FFV LSH
Source	in vitro E.coli expression system
Target Names	F46C5.2
Protein Names	Recommended name: Uncharacterized protein F46C5.2
Expression Region	20-202
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.