



Recombinant *Oryza sativa* subsp. japonica Copper transporter 4 (COPT4)

Product Code	CSB-CF606061OFG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q10KT6
Product Type	Transmembrane Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Sequence	MAMPMPMPPPGPGGDAPPAPTMMPGMAMPMTTGMSFTWGHRAVVLFP RW PGDRAGVGMYP LCLLLVLALAALAEALSAASRRLLDLLDLRSRGRRRRRRQQLLAAGVHA ARM GLAYLVM LAVMSFNAGVLLAAVAGHAAGFLLARSGLLGSRAAAPEIDGAAAAAAT SNGS SLHPSSE PKP
Source	in vitro E.coli expression system
Target Names	COPT4
Protein Names	Recommended name: Copper transporter 4 Short name= OsCOPT4
Expression Region	1-183
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.