



Recombinant African swine fever virus Envelope protein p54

Product Code	CSB-CF717644AEC
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q65272
Product Type	Transmembrane Protein
Immunogen Species	African swine fever virus (isolate Pig/Portugal/Lis 60/1960) (ASFV)
Sequence	MDSEFFQPVYPRHYGECLSPTSTPSFFSTHMYTILIAIVVLVIIIIVLIYLFSSRKKK AA AAIEEEDIQFINPYQDQQWAEVTPQPGTSKPAGATTASAGKPV TGRPATNRPA TNKPVTD NPVTDRLVMATGGPAAAPAAASAHPTEPYTTVTTQNTASQTMSAIENLRQRNT YTHKDLE NSL
Source	in vitro E.coli expression system
Protein Names	Recommended name: Envelope protein p54
Expression Region	1-183
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.