



Recombinant Danio rerio Mitochondrial inner membrane protease subunit 2 (immp2l)

Product Code	CSB-CF721486DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q6AZD4
Product Type	Transmembrane Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Sequence	MAQTGFGRRYFKAFVSGFFVAVPVTVTVLDRLAYVARVEGASMQPSLNPDGE SSPDVLL NRWSVRNYHVQRGDIVSVLSPKNPQQKIIKRVIIEGDFIKTLGYKNRYVRVPD GHLWIE GDHHGHSFDSNAFGPVSLGLVHGRASHIWPSPRWQRIEPSVPPDRRPLLNW DRAAEDKY DDD
Source	in vitro E.coli expression system
Target Names	immp2l
Protein Names	Recommended name: Mitochondrial inner membrane protease subunit 2 EC=3.4.21.- Alternative name(s): IMP2-like protein
Expression Region	1-183
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.