



Recombinant Human Fetal and adult testis-expressed transcript protein (FATE1)

Product Code	CSB-CF839274HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q969F0
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MAGGPPNTKAEMEMSLAEELNHGRQGENQEHLVIAEMMELGSRSRGASQKK QKLEQKAAG SASAKRVWNMTATRPKKMGSQPKPRMLRESGHGDAHLQEYAGNFQGIRFH YDRNPGTDA VAQTSLEEFNVLEMEVMRRQLYAVNRRLRALEEQGATWRHRETLIIAVLVSASI ANLWLW MNQ
Source	in vitro E.coli expression system
Target Names	FATE1
Protein Names	Recommended name: Fetal and adult testis-expressed transcript protein Alternative name(s): Cancer/testis antigen 43 Short name= CT43 Tumor antigen BJ-HCC-2
Expression Region	1-183
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.