



Recombinant Enterobacteria phage lambda Superinfection exclusion protein B (sieB)

Product Code	CSB-CF361146ECW
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P03762
Product Type	Transmembrane Protein
Immunogen Species	Escherichia phage lambda (Bacteriophage lambda)
Sequence	MMSIEMDPLVILGRVFSNEPLERTMYMIVIWVGLLLLSPDNWPEYVNERIGIPH VWHV FV FALAFSLAINVHRLSAIASARYKRFKLRKRIKMQNDKVRSVIQNLTEEQSMVLCA ALNEG RKYVVTSKQFPYISELIELGVLNKTFSRWNGKHILFPIEDIYWTELVASYDPYNIE IKPR PISK
Source	in vitro E.coli expression system
Target Names	sieB
Protein Names	Recommended name: Superinfection exclusion protein B
Expression Region	1-184
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.