



Recombinant Myeloproliferative leukemia virus Myeloproliferative leukemia protein (V-MPL)

Product Code	CSB-CF331243MMB
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P40931
Product Type	Transmembrane Protein
Immunogen Species	Myeloproliferative leukemia virus (MplV)
Sequence	LELRPRARYSLQLRARLNGPTYQGPWSAWSPPARVSTGSETAWITLVTALLLV LSLSALL GLLLLKWQFPAHYRRLRHALWPSLPDLHRVLGQYLRDTAALSPSKATVTDSCE EVEPSLL EILPKSSESTPLPLCPSQPQMDYRGLQPCLRTMPLSVCPPMAETGSCCTTHIA NHSYLPL SYWQ
Source	in vitro E.coli expression system
Target Names	V-MPL
Protein Names	Recommended name: Myeloproliferative leukemia protein
Expression Region	1-184
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.