



Recombinant Bovine Receptor expression-enhancing protein 6 (REEP6)

Product Code	CSB-CF653471BO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q32LG5
Product Type	Transmembrane Protein
Immunogen Species	Bos taurus (Bovine)
Sequence	MDGLRQRFERFLEQRNLATEALGALEAKTGVDKRYLATGAATLLSLYLLFGYG APLLCSL IGFAYPAYASIKAIESPskEDDTVWLTyWVYGLFGLAEFFSDLLLSWFPFYA GKCAFL LFCMAPGPWNGAHMLYHRIIRPLFLKHHEAVDSIVSDISGRALDVAAGMTKDA GKVSVNQ LQKAK
Source	in vitro E.coli expression system
Target Names	REEP6
Protein Names	Recommended name: Receptor expression-enhancing protein 6
Expression Region	1-185
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.