



Recombinant Mouse Receptor expression-enhancing protein 5 (Reep5)

Product Code	CSB-CF723351MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q60870
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MRERFDRFLHEKNCMTDLLAKLEAKTGVNRSFIALGVIGLVALYLVFGYGASLL CNLIGF GYPAYISMKAIESPNKDDDTQWLTYWVVYGVFSIAEFFSDLFLSWLPFYMLK CGFLLWC MAPSPANGAEMLYRRIIRPIFLRHESQVDSVVKDKAKETADAIKKEVKKAT VNLLGD VKKST
Source	in vitro E.coli expression system
Target Names	Reep5
Protein Names	Recommended name: Receptor expression-enhancing protein 5 Alternative name(s): GP106 Polyposis locus protein 1 homolog Protein TB2 homolog
Expression Region	1-185
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.