



Recombinant *Oryza sativa* subsp. *indica* Photosystem I assembly protein Ycf4 (ycf4)

Product Code	CSB-CF313588OFF
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0C514
Product Type	Transmembrane Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>indica</i> (Rice)
Sequence	MNWRSEHIWIELLKGSRKRGNFFWACILFLGSLGFLAVGASSYLGKNIISVLPS QQILFF PQGVVMSFYGIAGLFIAYLWCTILWNVGSGYDRFDRKEGVVCIFRWGFPGIK RRVFLRF LMRDIQSIRIQVKEGLFPRRILYMEIRGQGAIPLTRTDEKFFTPREIEQKAAELAY FLRI PMEVF
Source	in vitro E.coli expression system
Target Names	ycf4
Protein Names	Recommended name: Photosystem I assembly protein Ycf4
Expression Region	1-185
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.