



Recombinant Pig Selenoprotein S (SELS)

Product Code	CSB-CF515255PI
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	F1SR90
Product Type	Transmembrane Protein
Immunogen Species	Sus scrofa (Pig)
Sequence	MEQDGDQLSARPALETEGLRFLHVTVGSLLATYGWYIVFCCILLYVVFQKLSTR LRALRQ RHLDGAAAALPDVVVKRQEALAAARLKMQEELNAQVEKHKEKLRQLEEEKR RQKIERWD SVQEGRSYRGDARKRQEEDSPGPSTSSVIPKRKSDKKPLRGGGYNPLSGEG GGACSWRPG RRGPSGGUG
Source	in vitro E.coli expression system
Target Names	SELENOS
Protein Names	Recommended name: Selenoprotein S Short name= SelS
Expression Region	1-190
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.