



Recombinant Arabidopsis thaliana Oleosin 20.3 kDa (OL2)

Product Code	CSB-CF671459DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q42431
Product Type	Transmembrane Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Sequence	ANVDRDRRVHVDRTDKRVHQPNYEDDVGFGGYGGYGAGSDYKSRGPSTNQI LALIAGVPI GGTLLTAGLTLAGSVIGLLVSIPLFLLFSPVIVPAALTIGLAVTGILASGLFGLTGL SS VSWVLNYLRGTSDTVPEQLDYAKRRMADAVGYAGMKGKEMGQYVQDKAHE ARETEFMTET HEPGKARRGS
Source	in vitro E.coli expression system
Target Names	OL2
Protein Names	Recommended name: Oleosin 20.3 kDa Alternative name(s): Oleosin type 4
Expression Region	2-191
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.