



Recombinant Human Membrane protein FAM159A (FAM159A)

Product Code	CSB-CF744266HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q6UWV7
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MSGACTSYVSAEQEVVRGFSCPRPGGEAAAVFCCGFRDHKYCCDDPHSFFP YEHSYMWLWL SIGALIGLSVAAVLLAFIVTACVLCYLFISSKPHTKLDLGLSLQTAGPEEVSPDC QGVN TGMAAEVPKVSPLQQSYSCLN PQLESNEGQAVNSKRLHHCFMATVTTSDIP GSPEEASV PNPDL CGPVP
Source	in vitro E.coli expression system
Target Names	SHISAL2A
Protein Names	Recommended name: Membrane protein FAM159A
Expression Region	1-190
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.