



Recombinant Arabidopsis thaliana Protein PLANT CADMIUM RESISTANCE 10 (PCR10)

Product Code	CSB-CF851422DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8S8T8
Product Type	Transmembrane Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Sequence	MKEKKGHYVPPSYIPLTQSDADTEVETTPNLEIAVSESTKDDPRQWSSGICA CFDDMQS CCVGLFCPCYIFGKNAELLGSGTFAGPCLTHCISWALVNTICCFATNGALLGLP GCFVSC YACGYRKSLRAKYNLQEAPCGDFVTHFFCHLCAICQEYREIREQSSGSYPLDM KMAITNA PLAQTMESAN
Source	in vitro E.coli expression system
Target Names	PCR10
Protein Names	Recommended name: Protein PLANT CADMIUM RESISTANCE 10 Short name= AtPCR10
Expression Region	1-190
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.