



Recombinant Mouse Selenoprotein S (Sels)

Product Code	CSB-CF857718MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9BCZ4
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MDRDEEPLSARPALETESLRFLHVTVGSELLASYGWYILFSCILLYIVIQRSLRLR ALRQ RQLDQAETVLEPDVVVKRQEALAAARLRMQEDLNAQVEKHKEKLRQLEEEKR RQKIEMWD SMQEGRSYKRNSGRPQEEDGPGPSTSSVIPKGKSDKKPLRGGGYNPLTGEG GGTCSWRPG RRGPSGGUN
Source	in vitro E.coli expression system
Target Names	Selenos
Protein Names	Recommended name: Selenoprotein S Short name= SelS Alternative name(s): Minor histocompatibility antigen H47 VCP-interacting membrane protein
Expression Region	1-190
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.